


# Minimise your 6-8% blood sample re-collections

HealthRFID blood sample labelling solution gives state-of-the-art support to phlebotomists, nurses and ward staff to cut out error-prone, time-consuming and costly manual tasks. In addition, our solution avoids the inconvenience to the patient to have their blood sample retaken.

## Facts and Figures:

- **~40% of the time** for a traditional manual blood sample collection relates to confirming patient details and hand writing labels (especially when multiple tubes need to be labelled).
- Manual process steps often lead to up to **6-8% rejection and re-collection rates** because of incorrect, unreadable or missing label information.
- Manual processes rely on humans to follow guidelines. Our solution imposes best practice and compliance to Standard 5 as a built-in feature.
- Our solution ensures that blood sample collection and patient identification takes place at the bedside.
- We fully support 1D/2D barcodes and RFID technology.



HealthRFID blood sample labelling solution comes out of the box, can be used standalone and does not require any integration to other systems. Our solution provides intuitive workflows which can be fully customized and tailored to your existing procedures.

FOR MORE INFORMATION



+61 419 252 598



contact@healthrfid.com



HealthRFID

...the next generation of tracking